

## TLR2 Rabbit pAb Catalog Number: ES8213 For research use only.

## Overview

Product Name	TLR2 polyclonal antibody
Host species	Rabbit
Applications	ІНС-р
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/200-1/500
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous TLR2 protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	TLR2
Gene Name	TLR2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	7097
Human Swiss-Prot Number	O60603
Alternative Names	