



ELK Biotechnology

SOD1 Rabbit pAb
Catalog Number: ES8226
For research use only.

Overview

Product Name	SOD1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human SOD-1.
Specificity	SOD-1 Polyclonal Antibody detects endogenous levels of SOD-1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	SOD1
Gene Name	SOD1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	18kDa
Human Gene ID	6647
Human Swiss-Prot Number	P00441
Alternative Names	SOD1; Superoxide dismutase [Cu-Zn]; Superoxide dismutase 1; hSod1