



ELK Biotechnology

GAL4 DNA Binding Domain Rabbit pAb

Catalog Number: ES8231

For research use only.

Overview

Product Name	GAL4 DNA Binding Domain polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	
Recommended dilutions	Western Blot:1/10000-1/50000
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous GAL4 DNA Binding Domain protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GAL4 DNA Binding Domain
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	
Human Swiss-Prot Number	
Alternative Names	