



ELK Biotechnology

Histone H3 Rabbit pAb
Catalog Number: ES8234
For research use only.

Overview

Product Name	Histone H3 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/1000 Immunohistochemistry:1/200-1/500
Immunogen	Synthetic Peptide
Specificity	The antibody detects endogenous Histone H3 protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Histone H3
Gene Name	HIST1H3A; HIST1H3B; HIST1H3C; HIST1H3D; HIST1H3E; HIST1H3F; HIST1H3G; HIST1H3H; HIST1H3I; HIST1H3J
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	15-17kDa
Human Gene ID	8290
Human Swiss-Prot Number	P68431
Alternative Names	