



ELK Biotechnology

Stat2 Rabbit pAb
Catalog Number: ES8236
For research use only.

Overview

Product Name	Stat2 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/10000
Immunogen	Synthesized peptide derived from human Stat2 around the non-phosphorylation site of Y631.
Specificity	Stat2 Polyclonal Antibody detects endogenous levels of Stat2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Stat2
Gene Name	STAT2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	98kDa
Human Gene ID	6773
Human Swiss-Prot Number	P52630
Alternative Names	STAT2; Signal transducer and activator of transcription 2; p113