



ELK Biotechnology

MEK-5 Rabbit pAb
Catalog Number: ES8251
For research use only.

Overview

Product Name	MEK-5 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/1000-1/2000 Immunohistochemistry:1/200-1/500
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous MEK-5 protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	MEK-5
Gene Name	MAP2K5
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	49kDa
Human Gene ID	5067
Human Swiss-Prot Number	Q13163
Alternative Names	