



ELK Biotechnology

HAO1 Rabbit pAb
Catalog Number: ES8252
For research use only.

Overview

Product Name	HAO1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Mouse,Rat
Recommended dilutions	Western Blot:1/2000-1/5000
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous HAO1 protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	HAO1
Gene Name	HAO1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	41kDa
Human Gene ID	54363
Human Swiss-Prot Number	Q9UJM8
Alternative Names	