



ELK Biotechnology

Caspase-8 Rabbit pAb
Catalog Number: ES8254
For research use only.

Overview

Product Name	Caspase-8 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/1000-1/2000 Immunohistochemistry:1/50-1/100
Immunogen	Recombinant Protein
Specificity	The antibody detects endogenous Caspase-8 protein.
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Caspase-8
Gene Name	CASP8
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	43,57kDa
Human Gene ID	841
Human Swiss-Prot Number	Q14790
Alternative Names	