



# ELK Biotechnology

**AMPK $\beta$ 2 Rabbit pAb**  
**Catalog Number: ES8260**  
**For research use only.**

## Overview

---

<b>Product Name</b>	AMPK $\beta$ 2 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/1000-1/2000</b>
<b>Immunogen</b>	Recombinant Protein
<b>Specificity</b>	The antibody detects endogenous AMPK $\beta$ 2 protein.
<b>Formulation</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	AMPK $\beta$ 2
<b>Gene Name</b>	PRKAB2
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>30kDa</b>
<b>Human Gene ID</b>	5565
<b>Human Swiss-Prot Number</b>	O43741
<b>Alternative Names</b>	