



ELK Biotechnology

STAT5a Rabbit pAb
Catalog Number: ES8278
For research use only.

Overview

Product Name	STAT5a polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/20000
Immunogen	Synthesized peptide derived from human Stat5a around the non-phosphorylation site of S780.
Specificity	Stat5a Polyclonal Antibody detects endogenous levels of Stat5a protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	STAT5a
Gene Name	STAT5a
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	90kDa
Human Gene ID	6776
Human Swiss-Prot Number	P42229
Alternative Names	STAT5A; STAT5; Signal transducer and activator of transcription 5A