



# ELK Biotechnology

**TNF  $\alpha$  Rabbit pAb**  
**Catalog Number: ES8280**  
**For research use only.**

## Overview

---

Product Name	TNF $\alpha$ polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/50-1/200</b>
Immunogen	Recombinant Protein
Specificity	The antibody detects TNF $\alpha$
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	TNF $\alpha$
Gene Name	TNF $\alpha$
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>17,25kDa</b>
Human Gene ID	7124
Human Swiss-Prot Number	P01375
Alternative Names	