



# ELK Biotechnology

## eEF1A2 binding protein Rabbit pAb

Catalog Number: ES8281

For research use only.

### Overview

---

Product Name	eEF1A2 binding protein polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Rat,Mouse,Pig
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/50-1/200
Immunogen	Recombinant Protein
Specificity	The antibody detects eEF1A2 binding protein
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	eEF1A2 binding protein
Gene Name	eEF1A2 binding protein
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	115kDa
Human Gene ID	91156
Human Swiss-Prot Number	Q86VF2
Alternative Names	