



# ELK Biotechnology

**S100 Rabbit pAb**  
**Catalog Number: ES8282**  
**For research use only.**

## Overview

---

Product Name	S100 polyclonal antibody
Host species	Rabbit
Applications	<b>WB,IHC-p</b>
Species Cross-Reactivity	<b>Human,Mouse,Rat</b>
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>Immunohistochemistry:1/100-1/300</b> <b>ELISA:1/40000</b>
Immunogen	Recombinant Protein
Specificity	The antibody detects S100
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	S100
Gene Name	S100
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>15kDa</b>
Human Gene ID	6271
Human Swiss-Prot Number	P23297
Alternative Names	