



ELK Biotechnology

Estrogen Receptor α Rabbit pAb

Catalog Number: ES8285

For research use only.

Overview

Product Name	Estrogen Receptor α polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000
Immunogen	Recombinant Protein
Specificity	The antibody detects Estrogen Receptor α
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Estrogen Receptor α
Gene Name	Estrogen Receptor α
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	66kDa
Human Gene ID	2099
Human Swiss-Prot Number	P03372
Alternative Names	