



ELK Biotechnology

Flotillin-2 Rabbit pAb
Catalog Number: ES8286
For research use only.

Overview

Product Name	Flotillin-2 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Recombinant Protein
Specificity	The antibody detects Flotillin-2
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Flotillin-2
Gene Name	Flotillin-2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	49,58kDa
Human Gene ID	2319
Human Swiss-Prot Number	Q14254
Alternative Names	