



ELK Biotechnology

Bcl-XI Rabbit pAb
Catalog Number: ES8287
For research use only.

Overview

Product Name	Bcl-XI polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/50-1/200
Immunogen	Recombinant Protein
Specificity	The antibody detects Bcl-XI
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Bcl-XI
Gene Name	Bcl-XI
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	30kDa
Human Gene ID	598
Human Swiss-Prot Number	Q07817
Alternative Names	