



ELK Biotechnology

Nrf2 Rabbit pAb
Catalog Number: ES8288
For research use only.

Overview

| | |
|---------------------------------|--|
| Product Name | Nrf2 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,IHC-p |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/500-1/2000 ELISA:1/20000 |
| Immunogen | Synthesized peptide derived from human Nrf2 around the non-acetylation site of K599. |
| Specificity | Nrf2 Polyclonal Antibody detects endogenous levels of Nrf2 protein. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | Nrf2 |
| Gene Name | Nrf2 |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 75-100kDa |
| Human Gene ID | 4780 |
| Human Swiss-Prot Number | Q16236 |
| Alternative Names | NFE2L2; NRF2; Nuclear factor erythroid 2-related factor 2; NF-E2-related factor 2; NFE2-related factor 2; HEBP1; Nuclear factor, erythroid derived 2, like 2 |