



ELK Biotechnology

Bestrophin-1 Rabbit pAb
Catalog Number: ES8291
For research use only.

Overview

Product Name	Bestrophin-1 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthetic Peptide
Specificity	The antibody detects Bestrophin-1
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Bestrophin-1
Gene Name	Bestrophin-1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	67kDa
Human Gene ID	7439
Human Swiss-Prot Number	O76090
Alternative Names	