



ELK Biotechnology

Bombesin Receptor 2 Rabbit pAb

Catalog Number: ES8292

For research use only.

Overview

Product Name	Bombesin Receptor 2 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/40000
Immunogen	Synthetic Peptide
Specificity	The antibody detects Bombesin Receptor 2
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Bombesin Receptor 2
Gene Name	Bombesin Receptor 2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	50-70kDa
Human Gene ID	2925
Human Swiss-Prot Number	N/A
Alternative Names	