



ELK Biotechnology

CABP2 Rabbit pAb
Catalog Number: ES8295
For research use only.

Overview

| | |
|--------------------------|---|
| Product Name | CABP2 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,IHC-p |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/500-1/2000 Immunohistochemistry:1/100-1/300 ELISA:1/40000 |
| Immunogen | Synthetic Peptide |
| Specificity | The antibody detects CABP2 |
| Formulation | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | CABP2 |
| Gene Name | CABP2 |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 37kDa |
| Human Gene ID | 51475 |
| Human Swiss-Prot Number | Q9NPB3 |
| Alternative Names | |