



ELK Biotechnology

CaV α 2 δ 3 Rabbit pAb
Catalog Number: ES8303
For research use only.

Overview

Product Name	CaV α 2 δ 3 polyclonal antibody
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/1000-1/2000
Immunogen	Synthetic Peptide
Specificity	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CaV α 2 δ 3
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	120kDa
Human Gene ID	55799
Human Swiss-Prot Number	Q8I2S8
Alternative Names	