



# ELK Biotechnology

**CCKBR Rabbit pAb**  
**Catalog Number: ES8305**  
**For research use only.**

## Overview

---

<b>Product Name</b>	CCKBR polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/1000-1/2000</b>
<b>Immunogen</b>	Synthetic Peptide
<b>Specificity</b>	
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	CCKBR
<b>Gene Name</b>	
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>50kDa</b>
<b>Human Gene ID</b>	887
<b>Human Swiss-Prot Number</b>	P32239
<b>Alternative Names</b>	CCKBR; CCKRB; Gastrin/cholecystokinin type B receptor; CCK-B receptor; CCK-BR; Cholecystokinin-2 receptor; CCK2-R