



ELK Biotechnology

Connexin-26 Rabbit pAb
Catalog Number: ES8309
For research use only.

Overview

Product Name	Connexin-26 polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human, Mouse
Recommended dilutions	Immunohistochemistry: 1/100-1/200
Immunogen	Synthetic Peptide
Specificity	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Connexin-26
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	26kDa
Human Gene ID	2706
Human Swiss-Prot Number	P29033
Alternative Names	GJB2; Gap junction beta-2 protein; Connexin-26; Cx26