



# ELK Biotechnology

**EphA1 Rabbit pAb**  
**Catalog Number: ES8312**  
**For research use only.**

## Overview

---

<b>Product Name</b>	EphA1 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	WB,IHC-p
<b>Species Cross-Reactivity</b>	Human,Mouse,Rat
<b>Recommended dilutions</b>	<b>Western Blot:1/1000-1/2000</b> <b>Immunohistochemistry:1/100-1/200</b>
<b>Immunogen</b>	Synthetic Peptide
<b>Specificity</b>	
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	EphA1
<b>Gene Name</b>	
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>110kDa</b>
<b>Human Gene ID</b>	2041
<b>Human Swiss-Prot Number</b>	P21709
<b>Alternative Names</b>	EPHA1; EPH; EPHT; EPHT1; Ephrin type-A receptor 1; hEpha1; EPH tyrosine kinase; EPH tyrosine kinase 1; Erythropoietin-producing hepatoma receptor; Tyrosine-protein kinase receptor EPH