



# ELK Biotechnology

**KCNN2 (SK2) Rabbit pAb**  
**Catalog Number: ES8329**  
**For research use only.**

## Overview

---

Product Name	KCNN2 (SK2) polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	<b>Western Blot:1/1000-1/2000</b> <b>Immunohistochemistry:1/100-1/200</b>
Immunogen	Synthetic Peptide
Specificity	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	KCNN2 (SK2)
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>70,26kDa</b>
Human Gene ID	3781
Human Swiss-Prot Number	Q6PJI0
Alternative Names	