



ELK Biotechnology

KCNN3 (SK3) Rabbit pAb
Catalog Number: ES8330
For research use only.

Overview

Product Name	KCNN3 (SK3) polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Rat
Recommended dilutions	Immunohistochemistry:1/100-1/200
Immunogen	Synthetic Peptide
Specificity	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	KCNN3 (SK3)
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	82kDa
Human Gene ID	3782
Human Swiss-Prot Number	Q9UGI6
Alternative Names	