



ELK Biotechnology

Kv1.8 Rabbit pAb
Catalog Number: ES8339
For research use only.

Overview

| | |
|--------------------------|---|
| Product Name | Kv1.8 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,IHC-p |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/1000-1/2000 Immunohistochemistry:1/100-1/200 |
| Immunogen | Synthetic Peptide |
| Specificity | |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | Kv1.8 |
| Gene Name | |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 58kDa |
| Human Gene ID | 3744 |
| Human Swiss-Prot Number | Q16322 |
| Alternative Names | |