



ELK Biotechnology

Kv10.1 Rabbit pAb
Catalog Number: ES8340
For research use only.

Overview

Product Name	Kv10.1 polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/100-1/200
Immunogen	Synthetic Peptide
Specificity	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Kv10.1
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	50kDa
Human Gene ID	170850
Human Swiss-Prot Number	Q8TAE7
Alternative Names	KCNG3; Potassium voltage-gated channel subfamily G member 3; Voltage-gated potassium channel subunit Kv10.1; Voltage-gated potassium channel subunit Kv6.3