



# ELK Biotechnology

**DUSP6 Rabbit pAb**  
**Catalog Number: ES8344**  
**For research use only.**

## Overview

---

Product Name	DUSP6 polyclonal antibody
Host species	Rabbit
Applications	<b>WB</b>
Species Cross-Reactivity	<b>Human,Mouse,Rat</b>
Recommended dilutions	<b>Western Blot:1/1000-1/3000</b>
Immunogen	Synthetic Peptide
Specificity	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	DUSP6
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>42,44kDa</b>
Human Gene ID	1848
Human Swiss-Prot Number	Q16828
Alternative Names	DUSP6; MKP3; PYST1; Dual specificity protein phosphatase 6; Dual specificity protein phosphatase PYST1; Mitogen-activated protein kinase phosphatase 3; MAP kinase phosphatase 3; MKP-3