



ELK Biotechnology

TCEAL6 Rabbit pAb
Catalog Number: ES8359
For research use only.

Overview

Product Name	TCEAL6 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human TCEAL6. at AA range; 1-80
Specificity	TCEAL6 Polyclonal Antibody detects endogenous levels of TCEAL6 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	TCEAL6
Gene Name	TCEAL6
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	32kDa
Human Gene ID	85508
Human Swiss-Prot Number	Q6IPX3
Alternative Names	TCEAL6; Transcription elongation factor A protein-like 6; TCEA-like protein 6; Transcription elongation factor S-II protein-like 6