



ELK Biotechnology

TGF α Rabbit pAb
Catalog Number: ES8382
For research use only.

Overview

Product Name	TGF α polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human TGF α .
Specificity	TGF α Polyclonal Antibody detects endogenous levels of TGF α protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	TGF α
Gene Name	TGFA
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	15kDa
Human Gene ID	7039
Human Swiss-Prot Number	P01135
Alternative Names	TGFA; Protransforming growth factor alpha