



ELK Biotechnology

MAO-A Rabbit pAb
Catalog Number: ES8385
For research use only.

Overview

Product Name	MAO-A polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA:1/10000
Immunogen	Synthesized peptide derived from the Internal region of human MAO-A.
Specificity	MAO-A Polyclonal Antibody detects endogenous levels of MAO-A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	MAO-A
Gene Name	MAOA
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kDa
Human Gene ID	4128
Human Swiss-Prot Number	P21397
Alternative Names	MAOA; Amine oxidase [flavin-containing] A; Monoamine oxidase type A; MAO-A