



# ELK Biotechnology

**IL-4R $\alpha$  Rabbit pAb**  
**Catalog Number: ES8391**  
**For research use only.**

## Overview

---

Product Name	IL-4R $\alpha$ polyclonal antibody
Host species	Rabbit
Applications	<b>WB,ELISA</b>
Species Cross-Reactivity	<b>Human</b>
Recommended dilutions	<b>Western Blot:1/500-1/2000</b> <b>ELISA:1/10000</b>
Immunogen	Synthesized peptide derived from the Internal region of human IL-4R $\alpha$ .
Specificity	IL-4R $\alpha$ Polyclonal Antibody detects endogenous levels of IL-4R $\alpha$ protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	IL-4R $\alpha$
Gene Name	IL4R
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>90kDa</b>
Human Gene ID	3566
Human Swiss-Prot Number	P24394
Alternative Names	IL4R; IL4RA; 582J2.1; Interleukin-4 receptor subunit alpha; IL-4 receptor subunit alpha; IL-4R subunit alpha; IL-4R-alpha; IL-4RA; CD124