



ELK Biotechnology

β -actin Rabbit(A284) Rabbit pAb
Catalog Number: ES8450
For research use only.

Overview

Product Name	β -actin Rabbit(A284) polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/1000-1/2000 Immunohistochemistry:1/100-1/200
Immunogen	Synthetic Peptide of β -actin
Specificity	β -actin protein(A284) detects endogenous levels of ACTB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	β -actin Rabbit(A284)
Gene Name	ACTB
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	43kDa
Human Gene ID	60
Human Swiss-Prot Number	P60709
Alternative Names	ACTB