



ELK Biotechnology

PI3 Kinase P85 Rabbit Rabbit pAb

Catalog Number: ES8452

For research use only.

Overview

Product Name	PI3 Kinase P85 Rabbit polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/100-1/200
Immunogen	Recombinant Protein of PI3 Kinase P85 of PIK3R1
Specificity	PI3 Kinase P85 protein detects endogenous levels of PIK3R1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	PI3 Kinase P85 Rabbit
Gene Name	PIK3R1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	85kDa
Human Gene ID	5295
Human Swiss-Prot Number	P27986
Alternative Names	PIK3R1