



ELK Biotechnology

CD10 Rabbit Rabbit pAb
Catalog Number: ES8454
For research use only.

Overview

Product Name	CD10 Rabbit polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/100-1/200
Immunogen	Recombinant Protein of CD10 of MME
Specificity	CD10 protein detects endogenous levels of MME
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CD10 Rabbit
Gene Name	MME
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	85kDa
Human Gene ID	4311
Human Swiss-Prot Number	P08473
Alternative Names	MME