



ELK Biotechnology

NM23A Rabbit Rabbit pAb
Catalog Number: ES8455
For research use only.

Overview

Product Name	NM23A Rabbit polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/100-1/200
Immunogen	Recombinant Protein of NM23A of NME1
Specificity	NM23A protein detects endogenous levels of NME1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	NM23A Rabbit
Gene Name	NME1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	17kDa
Human Gene ID	4830
Human Swiss-Prot Number	P15531
Alternative Names	NME1