



# ELK Biotechnology

**ATM Rabbit Rabbit pAb**  
**Catalog Number: ES8456**  
**For research use only.**

## Overview

---

Product Name	ATM Rabbit polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/100-1/200
Immunogen	Recombinant Protein of ATM of ATM
Specificity	ATM protein detects endogenous levels of ATM
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ATM Rabbit
Gene Name	ATM
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	350kDa
Human Gene ID	472
Human Swiss-Prot Number	Q13315
Alternative Names	ATM