



ELK Biotechnology

MEK Kinase-1 Rabbit pAb
Catalog Number: ES8462
For research use only.

Overview

Product Name	MEK Kinase-1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from MEK Kinase-1 . at AA range; 1340-1420
Specificity	MEK Kinase-1 Polyclonal Antibody detects endogenous levels of MEK Kinase-1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	MEK Kinase-1
Gene Name	MAP3K1 MAPKKK1 MEKK MEKK1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	161kDa
Human Gene ID	4214
Human Swiss-Prot Number	Q13233
Alternative Names	mitogen-activated protein kinase kinase kinase 1