



ELK Biotechnology

TGF β RIII Rabbit pAb
Catalog Number: ES8473
For research use only.

Overview

Product Name	TGF β RIII polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from TGF β RIII at AA range; 731-780
Specificity	TGF β RIII Polyclonal Antibody detects endogenous levels of TGF β RIII
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	TGF β RIII
Gene Name	TGFBR3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	93kDa
Human Gene ID	7049
Human Swiss-Prot Number	Q03167
Alternative Names	transforming growth factor, beta receptor III