



ELK Biotechnology

Bcl-2 Rabbit pAb
Catalog Number: ES8474
For research use only.

Overview

Product Name	Bcl-2 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000 Immunohistochemistry:1/50-1/300
Immunogen	Synthesized peptide derived from Bcl-2 . at AA range; 10-90
Specificity	Bcl-2 Polyclonal Antibody detects endogenous levels of Bcl-2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Bcl-2
Gene Name	BCL2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	22kDa
Human Gene ID	596
Human Swiss-Prot Number	P10415
Alternative Names	B-cell CLL/lymphoma 2