



ELK Biotechnology

HLA Class I Rabbit pAb
Catalog Number: ES8483
For research use only.

Overview

Product Name	HLA Class I polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthetic Peptide of HLA Class I
Specificity	HLA Class I Polyclonal Antibody detects endogenous levels of HLA Class I
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	HLA Class I
Gene Name	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	42kDa
Human Gene ID	3105
Human Swiss-Prot Number	P01891/P01892/P04439/P05534/P10314/P10316/P13746/P16189/P18462/ P30443/P30450/P30453/P30455/P30456/P30457/P30512/Q09160
Alternative Names	