



ELK Biotechnology

Collagen IV Rabbit pAb
Catalog Number: ES8485
For research use only.

Overview

Product Name	Collagen IV polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from Collagen IV . at AA range; 1620-1670
Specificity	Collagen IV Polyclonal Antibody detects endogenous levels of Collagen IV
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Collagen IV
Gene Name	COL4A1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	130kDa
Human Gene ID	1282
Human Swiss-Prot Number	P02462
Alternative Names	collagen, type IV, alpha 1