



ELK Biotechnology

NF- κ B p65 Rabbit pAb
Catalog Number: ES8488
For research use only.

Overview

Product Name	NF- κ B p65 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthetic Peptide of NF- κ B p65
Specificity	NF- κ B p65 Polyclonal Antibody detects endogenous levels of NF- κ B p65
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	NF- κ B p65
Gene Name	RELA NFKB3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	60kDa
Human Gene ID	5970
Human Swiss-Prot Number	Q04206
Alternative Names	v-rel reticuloendotheliosis viral oncogene homolog A (avian)