



ELK Biotechnology

p53 Rabbit pAb
Catalog Number: ES8490
For research use only.

Overview

Product Name	p53 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthetic Peptide of p53
Specificity	p53 Polyclonal Antibody detects endogenous levels of p53
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	p53
Gene Name	TP53 P53
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	53kDa
Human Gene ID	7157
Human Swiss-Prot Number	P04637
Alternative Names	tumor protein p53