



FAM3D rabbit pAb

Cat No.:ES8496

For research use only


Overview

Product Name	FAM3D rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000, ELISA 1:10000-20000
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human FAM3D. AA range:121-170
Specificity	FAM3D Polyclonal Antibody detects endogenous levels of FAM3D
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	FAM3D
Gene Name	FAM3D
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	22kD
Human Gene ID	131177
Human Swiss-Prot Number	Q96BQ1
Alternative Names	Protein FAM3D
Background	similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small intestine.,





ELK Biotechnology



Western Blot analysis of HeLa cells using FAM3D Polyclonal Antibody diluted at 1:500. Secondary antibody(catalog#:RS0002) was diluted at 1:20000



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C