



ELK Biotechnology

Casein Kinase I α Rabbit pAb
Catalog Number: ES8509
For research use only.

Overview

Product Name	Casein Kinase I α polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from Casein Kinase I α at AA range; 161-210
Specificity	Casein Kinase I α Polyclonal Antibody detects endogenous levels of Casein Kinase I α
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Casein Kinase I α
Gene Name	CSNK1A1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	40kDa
Human Gene ID	1452
Human Swiss-Prot Number	P48729
Alternative Names	Casein kinase I isoform alpha (CKI-alpha) (EC 2.7.11.1) (CK1)