



ELK Biotechnology

Caldesmon Rabbit pAb
Catalog Number: ES8513
For research use only.

Overview

Product Name	Caldesmon polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from Caldesmon at AA range; 1-50
Specificity	Caldesmon Polyclonal Antibody detects endogenous levels of Caldesmon
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Caldesmon
Gene Name	CALD1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	95kDa
Human Gene ID	800
Human Swiss-Prot Number	Q05682
Alternative Names	Caldesmon (CDM)