



ELK Biotechnology

Vimentin Rabbit pAb
Catalog Number: ES8514
For research use only.

Overview

Product Name	Vimentin polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from Vimentin at AA range; 371-420
Specificity	Vimentin Polyclonal Antibody detects endogenous levels of Vimentin
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Vimentin
Gene Name	VIM
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	55kDa
Human Gene ID	7431
Human Swiss-Prot Number	P08670
Alternative Names	Vimentin