



ELK Biotechnology

α -SMA Rabbit pAb
Catalog Number: ES8515
For research use only.

Overview

Product Name	α -SMA polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000 Immunohistochemistry:1/50-1/300
Immunogen	Synthesized peptide derived from α -SMA at AA range; 171-220
Specificity	α -SMA Polyclonal Antibody detects endogenous levels of α -SMA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	α -SMA
Gene Name	ACTA2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	42kDa
Human Gene ID	
Human Swiss-Prot Number	P62736/P68133/P63267/P63261/P68032
Alternative Names	